



# RAPGEF2 Antibody, Biotin conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA897493LD01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q9Y4G8   |
| <b>Immunogen</b>           | Recombinant Human Rap guanine nucleotide exchange factor 2 protein (1393-1498AA)   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Biotin   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4   |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Rap guanine nucleotide exchange factor 2 (Cyclic nucleotide ras GEF) (CNrasGEF) (Neural RAP guanine nucleotide exchange protein) (nRap GEP) (PDZ domain-containing guanine nucleotide exchange factor 1) (PDZ-GEF1) (RA-GEF-1) (Ras/Rap1-associating GEF-1), RAPGEF2, KIAA0313 NRAPGEP PDZGEF1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Neuroscience   |
| <b>Target Names</b>        | RAPGEF2  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.  |