




Butyryl-HIST1H3A (K9) Antibody

Product Code	CSB-PA010418OA09butHU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431
Immunogen	Peptide sequence around site of Butyryl-Lys (9) derived from Human Histone H3.1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIST1H3A

Image



Western Blot

Positive WB detected in: Hela whole cell lysate(30µg, treated with 30mM sodium butyrate for 4h), HepG2 whole cell lysate(30µg,treated with 30mM sodium butyrate for 4h), Hela whole cell lysate(30µg),HepG2 whole cell lysate(30µg)
 All lanes: HIST1H3A antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/20000 dilution
 Predicted band size: 15-18 kDa
 Observed band size: 16 kDa
 Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.