




Recombinant Human RNA-binding protein RO60 (RO60), partial

Product Code	CSB-EP025071HU1
Relevance	RNA-binding protein that binds to misfolded non-coding RNAs, pre-5S rRNA, and several small Cytoplasmic domain RNA molecules known as Y RNAs. May stabilize some of these RNAs and protect th from degradation.
Abbreviation	Recombinant Human RO60 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10155
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ESVNQMQLNEKQIANSQDGYVWQVTDNMNRLHRFLCFGSEGGTYIKEQKL GLENAEALIRLIEDGRGCEVIQEIKSFSQEGRTTKQEPMLFALAICSQCSDISTK QAAFKAVSEVCRIPHTLFTFIQFKKDLKESMKCGMWGRALRKAADWYNEKGG MALALAVTKYKQRNGWSHKDLLRLSHLKPSSSEGLAIVTKYITKGWKEVHELYK EKALSVETEKLLKYLEAVEKVKRTRDELEVIHLIEEHRLVREHLLTNHLKSKEVW KALLQEMPLTALLRNLGKMTANSVLEPGNSEVSLVCEKLCNEKLLKKARIHPFH ILIALETYKTGHGLRGKLKWRPDEEILKALDAAFYKTFKTVEPTGKRFLLAVDVS ASMNQRVLGSILNASTVAAAMCMVTRTEKDSYVVAFSDEMVPVPCPVTTDRTL QQVLMAMSQIPAGGTDCSLPMIWAQKTNTPADVFIVFTDNETFAGGVHPIAL REYRKKMDIPAKLIVCGMTSNGFTIADPDDRGMLDMCGFDTGALDVIRNFTL
Source	E.coli
Target Names	TROVE2
Expression Region	3-535aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	64.1kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.