



Recombinant Human Collagen alpha-1 (XIII) chain (COL13A1), partial

Product Code	CSB-BP735980HU1
Abbreviation	COL13A1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TAT6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	KGSKGEPGKGMVDYNGNINEALQEIRTLALMGPPGLPGQIGPPGAPGIPGQK GEIGLPGPPGHDGEKGRGKPGDMGPPGPQGGPKDGPPGVKGENGHPGS PGEKGEKGETGQAGSPGEKGEAGEKGNPGAIEVPLPGPEGPPGPPGLQGV PGPKGEAGLDGAKGEKGFQGEKDRGPLGLPGASGLDGRPGPPGTPGPIGV PGPAGPKGERGSKGDPGMTGPTGAAGLPLHGPPGDKGNRGERGKKGSRG PKGDKGDQGAPGLDAPCPLGEDGLPVQGCWNK
Source	Baculovirus
Target Names	COL13A1
Protein Names	Recommended name: Collagen alpha-1(XIII) chain Alternative name(s): COLXIII A1
Expression Region	432-717aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.