



SPTA1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA958244 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02549 |
| Immunogen | Synthetic peptide of Human SPTA1 |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:60-1:450 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.05% sodium azide, 40% glycerol, pH7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SPTA1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |