



SPTA1 Antibody

Product Code	CSB-PA966552
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02549
Immunogen	Synthetic peptide of Human SPTA1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.05% sodium azide, 40% glycerol, pH7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPTA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.