




Recombinant Bordetella bronchiseptica Azurin

Product Code	CSB-YP358313BTYb1
Abbreviation	Recombinant Bordetella bronchiseptica Azurin protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A321
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Bordetella bronchiseptica (strain ATCC BAA-588 / NCTC 13252 / RB50) (Alcaligenes bronchisepticus)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	AECSVDIAGTDQMQFDKKAIEVSKSCKQFTVNLKHTGKLPRNVMGHNWVLTK TADMQAVEKDGIAAGLDNQYLKAGDTRVLAHTKVLGGGESDSVTFDVAKLAA GDDYTFFCSFPGHGALMKGTCLKLVD
Research Area	Others
Source	Yeast
Target Names	BB3856
Expression Region	22-150aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	17.8 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.