



CLCA1 Antibody, Biotin conjugated

Product Code	CSB-PA005474LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8K7I4
Immunogen	Recombinant Human Calcium-activated chloride channel regulator 1 protein (715-785AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-activated chloride channel regulator 1 (EC 3.4.-.-) (Calcium-activated chloride channel family member 1) (hCLCA1) (Calcium-activated chloride channel protein 1) (CaCC-1) (hCaCC-1), CLCA1, CACC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CLCA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.