



CYP7A1 Antibody, Biotin conjugated

Product Code	CSB-PA006460LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22680
Immunogen	Recombinant Human Cholesterol 7-alpha-monooxygenase protein (303-423AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cholesterol 7-alpha-monooxygenase (EC 1.14.14.23) (CYPVII) (Cholesterol 7-alpha-hydroxylase) (Cytochrome P450 7A1), CYP7A1, CYP7
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP7A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.