



TCAF2 Antibody, FITC conjugated

Product Code	CSB-PA008015LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NFQ2
Immunogen	Recombinant Human TRPM8 channel-associated factor 2 protein (472-590AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TRPM8 channel-associated factor 2 (TRP channel-associated factor 2), TCAF2, FAM115C FAM139A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TCAF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.