



TCAF2 Antibody, Biotin conjugated

Product Code	CSB-PA008015LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NFQ2
Immunogen	Recombinant Human TRPM8 channel-associated factor 2 protein (472-590AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TRPM8 channel-associated factor 2 (TRP channel-associated factor 2), TCAF2, FAM115C FAM139A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TCAF2