




FAM13A Antibody

Product Code	CSB-PA008075LA01HU
Abbreviation	Protein FAM13A
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94988
Immunogen	Recombinant Human Protein FAM13A protein (1-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:200-1:500
Relevance	cytosol, GTPase activator activity, regulation of small GTPase mediated signal transduction
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM13A, FAM13A, FAM13A1 KIAA0914
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAM13A

Image



IHC image of CSB-PA008075LA01HU diluted at 1:300 and staining in paraffin-embedded human prostate cancer performed on a Leica Bond™ system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.