



FAM13A Antibody, HRP conjugated

Product Code	CSB-PA008075LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94988
Immunogen	Recombinant Human Protein FAM13A protein (1-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM13A, FAM13A, FAM13A1 KIAA0914
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAM13A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.