



FEZF1 Antibody, Biotin conjugated

Product Code	CSB-PA008602LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0PJY2
Immunogen	Recombinant Human Fez family zinc finger protein 1 protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fez family zinc finger protein 1 (Zinc finger protein 312B), FEZF1, FEZ ZNF312B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FEZF1