



FEZF1 Antibody, Biotin conjugated

Product Code	CSB-PA008602LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0PJY2
Immunogen	Recombinant Human Fez family zinc finger protein 1 protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fez family zinc finger protein 1 (Zinc finger protein 312B), FEZF1, FEZ ZNF312B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FEZF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.