



GCLC Antibody, HRP conjugated

Product Code	CSB-PA009321LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48506
Immunogen	Recombinant Human Glutamate--cysteine ligase catalytic subunit protein (78-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate--cysteine ligase catalytic subunit (EC 6.3.2.2) (GCS heavy chain) (Gamma-ECS) (Gamma-glutamylcysteine synthetase), GCLC, GLCL GLCLC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GCLC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.