



GGCT Antibody, HRP conjugated

Product Code	CSB-PA009387LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75223
Immunogen	Recombinant Human Gamma-glutamylcyclotransferase protein (54-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-glutamylcyclotransferase (EC 4.3.2.9) (Cytochrome c-releasing factor 21), GGCT, C7orf24 CRF21
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GGCT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.