



GGCT Antibody, Biotin conjugated

Product Code	CSB-PA009387LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75223
Immunogen	Recombinant Human Gamma-glutamylcyclotransferase protein (54-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-glutamylcyclotransferase (EC 4.3.2.9) (Cytochrome c-releasing factor 21), GGCT, C7orf24 CRF21
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GGCT