



GPN3 Antibody, Biotin conjugated

Product Code	CSB-PA009726LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHW5
Immunogen	Recombinant Human GPN-loop GTPase 3 protein (181-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPN-loop GTPase 3 (ATP-binding domain 1 family member C), GPN3, ATPBD1C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.