



HMGCS1 Antibody, FITC conjugated

Product Code	CSB-PA010566LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01581
Immunogen	Recombinant Human Hydroxymethylglutaryl-CoA synthase, cytoplasmic protein (360-520AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxymethylglutaryl-CoA synthase, cytoplasmic (HMG-CoA synthase) (EC 2.3.3.10) (3-hydroxy-3-methylglutaryl coenzyme A synthase), HMGCS1, HMGCS
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HMGCS1