



KIR2DL2 Antibody, HRP conjugated

Product Code	CSB-PA012353LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43627
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DL2 protein (22-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 2DL2 (CD158 antigen-like family member B1) (MHC class I NK cell receptor) (Natural killer-associated transcript 6) (NKAT-6) (p58 natural killer cell receptor clone CL-43) (p58 NK receptor CL-43) (CD antigen CD158b1), KIR2DL2, CD158B1 NKAT6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR2DL2