



KIR3DL2 Antibody, FITC conjugated

Product Code	CSB-PA012365LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43630
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 3DL2 protein (361-455AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 3DL2 (CD158 antigen-like family member K) (MHC class I NK cell receptor) (Natural killer-associated transcript 4) (NKAT-4) (p70 natural killer cell receptor clone CL-5) (p70 NK receptor CL-5) (CD antigen CD158k), KIR3DL2, CD158K NKAT4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR3DL2