



Krt18 Antibody, HRP conjugated

Product Code	CSB-PA012519LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05784
Immunogen	Recombinant Mouse Keratin, type I cytoskeletal 18 protein (200-423AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 18 (Cytokeratin endo B) (Keratin D) (Cytokeratin-18) (CK-18) (Keratin-18) (K18), Krt18, Kerd Krt1-18
Immunogen Species	Mus musculus (Mouse)
Research Area	Tags & Cell Markers
Target Names	Krt18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.