



KRT76 Antibody, Biotin conjugated

Product Code	CSB-PA012570LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01546
Immunogen	Recombinant Human Keratin, type II cytoskeletal 2 oral protein (173-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 2 oral (Cytokeratin-2P) (CK-2P) (K2P) (Keratin-76) (K76) (Type-II keratin Kb9), KRT76, KRT2B KRT2P
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	KRT76
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.