




LEP Antibody

Product Code	CSB-PA012870LA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50595
Immunogen	Recombinant Bovine Leptin protein (22-167AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leptin (Obesity factor), LEP, OB OBS
Immunogen Species	Bos taurus (Bovine)
Research Area	Signal Transduction
Target Names	LEP

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: LEP antibody at 4.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 21 kDa

Observed band size: 21 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.