



LEP Antibody, FITC conjugated

Product Code	CSB-PA012870LC01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50595
Immunogen	Recombinant Bovine Leptin protein (22-167AA)
Raised In	Rabbit
Species Reactivity	Bovine
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leptin (Obesity factor), LEP, OB OBS
Immunogen Species	Bos taurus (Bovine)
Research Area	Signal Transduction
Target Names	LEP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.