



LEP Antibody, FITC conjugated

Product Code	CSB-PA012870LC01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50595
Immunogen	Recombinant Bovine Leptin protein (22-167AA)
Raised In	Rabbit
Species Reactivity	Bovine
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leptin (Obesity factor), LEP, OB OBS
Immunogen Species	Bos taurus (Bovine)
Research Area	Signal Transduction
Target Names	LEP