



LGR5 Antibody, Biotin conjugated

Product Code	CSB-PA012906LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75473
Immunogen	Recombinant Human Leucine-rich repeat-containing G-protein coupled receptor 5 protein (241-352AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing G-protein coupled receptor 5 (G-protein coupled receptor 49) (G-protein coupled receptor 67) (G-protein coupled receptor HG38), LGR5, GPR49 GPR67
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LGR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.