



LGR5 Antibody, Biotin conjugated

Product Code	CSB-PA012906LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75473
Immunogen	Recombinant Human Leucine-rich repeat-containing G-protein coupled receptor 5 protein (241-352AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing G-protein coupled receptor 5 (G-protein coupled receptor 49) (G-protein coupled receptor 67) (G-protein coupled receptor HG38), LGR5, GPR49 GPR67
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LGR5