



MDH1 Antibody, Biotin conjugated

Product Code	CSB-PA013621ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40925
Immunogen	Recombinant Human Malate dehydrogenase, cytoplasmic protein (1-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malate dehydrogenase, cytoplasmic (EC 1.1.1.37) (Cytosolic malate dehydrogenase) (Diiodophenylpyruvate reductase) (EC 1.1.1.96), MDH1, MDHA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MDH1