




MEX3A Antibody

Product Code	CSB-PA013740LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1L020
Immunogen	Recombinant Human RNA-binding protein MEX3A protein (416-466AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein MEX3A (RING finger and KH domain-containing protein 4), MEX3A, RKHD4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MEX3A

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, SH-SY5Y whole cell lysate


All lanes: MEX3A antibody at 4.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa



IHC image of CSB-PA013740LA01HU diluted at 1:300 and staining in paraffin-embedded human brain tissue performed on a Leica BondTM system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.