



MEX3A Antibody, Biotin conjugated

Product Code	CSB-PA013740LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1L020
Immunogen	Recombinant Human RNA-binding protein MEX3A protein (416-466AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein MEX3A (RING finger and KH domain-containing protein 4), MEX3A, RKHD4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MEX3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.