



PHF7 Antibody, Biotin conjugated

Product Code	CSB-PA017918LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWX1
Immunogen	Recombinant Human PHD finger protein 7 protein (277-368AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PHD finger protein 7 (Testis development protein NYD-SP6), PHF7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PHF7