



PIGN Antibody, HRP conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA017977LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O95427 |
| Immunogen | Recombinant Human GPI ethanolamine phosphate transferase 1 protein (27-181AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | GPI ethanolamine phosphate transferase 1 (EC 2.-.-.-) (MCD4 homolog) (Phosphatidylinositol-glycan biosynthesis class N protein) (PIG-N), PIGN, MCD4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | PIGN |