



PIGN Antibody, HRP conjugated

Product Code	CSB-PA017977LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95427
Immunogen	Recombinant Human GPI ethanolamine phosphate transferase 1 protein (27-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPI ethanolamine phosphate transferase 1 (EC 2.-.-.) (MCD4 homolog) (Phosphatidylinositol-glycan biosynthesis class N protein) (PIG-N), PIGN, MCD4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PIGN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.