



Prl Antibody, FITC conjugated

Product Code	CSB-PA018724LC01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01237
Immunogen	Recombinant Rat Prolactin protein (30-226AA)
Raised In	Rabbit
Species Reactivity	Rat
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolactin (PRL), Prl
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Prl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.