



Prl Antibody, Biotin conjugated

Product Code	CSB-PA018724LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01237
Immunogen	Recombinant Rat Prolactin protein (30-226AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolactin (PRL), Prl
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Prl