



RBP3 Antibody

Product Code	CSB-PA019482LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10745
Immunogen	Recombinant Human Retinol-binding protein 3 protein (1146-1246AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinol-binding protein 3 (Interphotoreceptor retinoid-binding protein) (IRBP) (Interstitial retinol-binding protein), RBP3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RBP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.