



RYR3 Antibody, Biotin conjugated

Product Code	CSB-PA020621LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15413
Immunogen	Recombinant Human Ryanodine receptor 3 protein (987-1147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ryanodine receptor 3 (RYR-3) (RyR3) (Brain ryanodine receptor-calcium release channel) (Brain-type ryanodine receptor) (Type 3 ryanodine receptor), RYR3, HBRR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RYR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.