



RYR3 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA020621LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q15413 |
| Immunogen | Recombinant Human Ryanodine receptor 3 protein (987-1147AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ryanodine receptor 3 (RYR-3) (RyR3) (Brain ryanodine receptor-calcium release channel) (Brain-type ryanodine receptor) (Type 3 ryanodine receptor), RYR3, HBRR |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | RYR3 |