



SERPINB2 Antibody, HRP conjugated

Product Code	CSB-PA021070LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05120
Immunogen	Recombinant Human Plasminogen activator inhibitor 2 protein (293-389AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasminogen activator inhibitor 2 (PAI-2) (Monocyte Arg-serpin) (Placental plasminogen activator inhibitor) (Serpine B2) (Urokinase inhibitor), SERPINB2, PAI2 PLANH2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SERPINB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.