



SERPINB2 Antibody, Biotin conjugated

Product Code	CSB-PA021070LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05120
Immunogen	Recombinant Human Plasminogen activator inhibitor 2 protein (293-389AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasminogen activator inhibitor 2 (PAI-2) (Monocyte Arg-serpin) (Placental plasminogen activator inhibitor) (Serpine B2) (Urokinase inhibitor), SERPINB2, PAI2 PLANH2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SERPINB2