



SH3BP2 Antibody, Biotin conjugated

Product Code	CSB-PA021224LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78314
Immunogen	Recombinant Human SH3 domain-binding protein 2 protein (165-301AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-binding protein 2 (3BP-2), SH3BP2, 3BP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SH3BP2