



SLC1A7 Antibody, Biotin conjugated

Product Code	CSB-PA021438ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00341
Immunogen	Recombinant Human Excitatory amino acid transporter 5 protein (115-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Excitatory amino acid transporter 5 (Retinal glutamate transporter) (Solute carrier family 1 member 7), SLC1A7, EAAT5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC1A7