



SLC35D3 Antibody, FITC conjugated

Product Code	CSB-PA021596LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5M8T2
Immunogen	Recombinant Human Solute carrier family 35 member D3 protein (304-416AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 35 member D3 (Fringe connection-like protein 1), SLC35D3, FRCL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC35D3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.