



SRGN Antibody, FITC conjugated

Product Code	CSB-PA022664EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10124
Immunogen	Recombinant Human Serglycin protein (28-154AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serglycin (Hematopoietic proteoglycan core protein) (Platelet proteoglycan core protein) (P.PG) (Secretory granule proteoglycan core protein), SRGN, PRG PRG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SRGN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.