



SULT1A2 Antibody, FITC conjugated

Product Code	CSB-PA022934LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50226
Immunogen	Recombinant Human Sulfotransferase 1A2 protein (195-284AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sulfotransferase 1A2 (ST1A2) (EC 2.8.2.1) (Aryl sulfotransferase 2) (Phenol sulfotransferase 2) (Phenol-sulfating phenol sulfotransferase 2) (P-PST 2), SULT1A2, STP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SULT1A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.