



WIZ Antibody, Biotin conjugated

Product Code	CSB-PA026123LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95785
Immunogen	Recombinant Human Protein Wiz protein (1469-1574AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Wiz (Widely-interspaced zinc finger-containing protein) (Zinc finger protein 803), WIZ, ZNF803
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	WIZ