



TAP2 Antibody, HRP conjugated

Product Code	CSB-PA22539B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03519
Immunogen	Recombinant Human Antigen peptide transporter 2 protein (473-615AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Antigen peptide transporter 2 (APT2) (ATP-binding cassette sub-family B member 3) (Peptide supply factor 2) (Peptide transporter PSF2) (PSF-2) (Peptide transporter TAP2) (Peptide transporter involved in antigen processing 2) (Really interesting new gene 11 protein), TAP2, ABCB3 PSF2 RING11 Y1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.