



DNAJC15 Antibody, Biotin conjugated

Product Code	CSB-PA22849D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5T4
Immunogen	Recombinant Human DnaJ homolog subfamily C member 15 protein (59-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 15 (Cell growth-inhibiting gene 22 protein) (Methylation-controlled J protein) (MCJ), DNAJC15, DNAJD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJC15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.